

# Quick Reference Guide

## ViewRNA™ ISH Tissue 2-Plex Assay

Fresh Frozen or OCT-Embedded Frozen Samples

**!** **IMPORTANT:** If running ViewRNA ISH Tissue 2-Plex Assay for the first time, please refer to the *ViewRNA ISH Tissue 2-Plex Assay User Manual* to review assay guidelines and detailed procedures.

### Part 1: Sample Preparation and Target Probe Hybridization

Step	Task
1 Prepare and Chill 10% NBF	<ul style="list-style-type: none"><li>■ Prepare 200 mL of 10% NBF in 1X PBS and chill at 4 °C for 1 hr</li></ul>
2 Fix Tissue Overnight 16-18 hr	<ul style="list-style-type: none"><li>■ Transfer frozen tissue slides into chilled 10% NBF and incubate at 4 °C for 16-18 hr</li></ul>
3 Preparation for Part 1	<ul style="list-style-type: none"><li>■ Verify hybridization system is set to 40 °C and appropriately humidified</li><li>■ Prepare:<ul style="list-style-type: none"><li>□ 2 L 1X PBS</li><li>□ 200 mL 50% ethanol</li><li>□ 200 mL 70% ethanol</li><li>□ 4 L Wash Buffer</li><li>□ 200 mL Storage Buffer</li></ul></li><li>■ Ensure availability of:<ul style="list-style-type: none"><li>□ 1 L ddH<sub>2</sub>O</li><li>□ 200 mL Gill's Hematoxylin I</li><li>□ 200 mL of 3 µg/mL DAPI in 1X PBS (optional, for fluorescence detection)</li></ul></li><li>■ Thaw probe set(s) on ice</li><li>■ Prewarm 40 mL 1X PBS and Probe Set Diluent QF to 40 °C</li><li>■ Optional for 1-day assay:<ul style="list-style-type: none"><li>□ Prewarm PreAmplifier Mix QT, Amplifier Mix QT, and Label Probe Diluent QF to 40 °C</li><li>□ Briefly spin down Label Probe 1-AP, Label probe 6-AP, and Blue Reagents, place on ice</li><li>□ Bring Fast Red Tablets, Napthol Buffer, and AP Enhancer Solution to RT</li><li>□ Prepare 1 L 0.01% ammonium hydroxide under a fume hood</li></ul></li></ul>
4 Wash Slides 5 min	Wash slides 2 times in 1X PBS, 1 min each wash
5 Tissue Dehydration 30 min	<ul style="list-style-type: none"><li>■ Sequentially soak the rack of slides in 50%, 70%, and then 100% ethanol at RT without agitation, 10 min each time</li><li>■ Remove slides from 100% ethanol and drain on a paper towel.</li><li>■ Bake entire slide rack at 60 °C for 60 min</li></ul>
6 Draw Hydrophobic Barrier 40 min	<ul style="list-style-type: none"><li>■ Create hydrophobic barrier</li><li>■ Allow slides to air dry at RT for 20-30 min</li></ul>
7 Protease Digestion and Fixation 30-50 min, depending on optimized time	<ul style="list-style-type: none"><li>■ Prepare 1:100 working protease solution in prewarmed 1X PBS</li><li>■ Add working protease solution to slides</li><li>■ Incubate at 40 °C for optimal time determined in the pretreatment assay optimization procedure</li><li>■ Wash slides 2 times in 1X PBS, 1 min each wash</li><li>■ Fix slides in 10% NBF at RT for 5 min under a fume hood</li><li>■ Wash slides 2 times in 1X PBS, 1 min each wash</li></ul>
8 Target Probe Set Hybridization 2 hr 10 min	<ul style="list-style-type: none"><li>■ Prepare 1:40 working probe set solution in prewarmed Probe Set Diluent QT</li><li>■ Add working probe set solution to slides</li><li>■ Incubate at 40 °C for 2 hr</li></ul>
9 Wash Slides 8 min	Wash slides 3 times in Wash Buffer, 2 min each wash

Step	Task
<b>10 Optional Stop Point</b> 1 min	<ul style="list-style-type: none"> <li>■ Store slides in Storage Buffer at RT for up to 24 hr. Cover dish with lid or sealing film to prevent evaporation.</li> <li>■ Store 1X PBS and Wash Buffer at RT for use in Part 2</li> </ul>

## Part 2: Signal Amplification and Detection

Step	Task
<b>11 Preparation for Part 2</b> 10 min	<ul style="list-style-type: none"> <li>■ Pour Gill's Hematoxylin I into a clear staining dish, store at RT protected from light.</li> <li>■ Optional – Prepare 200 mL 3 µg/mL DAPI, store at 4 °C until use or place on ice</li> <li>■ Prewarm PreAmplifier Mix QT, Amplifier Mix QT and Label Probe Diluent QF to 40 °C</li> <li>■ Place Label Probe 1-AP, Label Probe 6-AP, and Blue reagents on ice</li> <li>■ Bring Fast Red Tablets, Naphthol Buffer, AP Enhancer, and Blue Buffer to RT</li> <li>■ Prepare 1 L 0.01% ammonium hydroxide in ddH<sub>2</sub>O</li> </ul>
<b>12 Wash Slides</b> 8 min	Remove slides from Storage Buffer and wash 3 times with Wash Buffer with constant and vigorous agitation, 2 min each wash
<b>13 PreAmplifier Hybridization</b> 35 min	Add PreAmplifier Mix QT directly to slides, incubate at 40 °C for 25 min
<b>14 Wash Slides</b> 8 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 2 min each wash
<b>15 Amplifier Hybridization</b> 20 min	Add Amplifier Mix QT directly to slides, incubate at 40 °C for 15 min
<b>16 Wash Slides</b> 8 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 2 min each wash
<b>17 Label Probe 6-AP Hybridization</b> 20 min	<ul style="list-style-type: none"> <li>■ Prepare 1:1000 working Label Probe 6-AP solution</li> <li>■ Add working Label Probe 6-AP solution to slides, incubate at 40 °C for 15 min</li> </ul>
<b>18 Wash Slides</b> 12 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 3 min each wash
<b>19 Apply Fast Blue Substrate</b> 35 min	<ul style="list-style-type: none"> <li>■ Prepare Fast Blue Substrate</li> <li>■ Add Fast Blue Substrate to slides, incubate at RT in the dark for 30 min</li> </ul>
<b>21 Wash Slides</b> 12 min	Wash slides 3 times in Wash Buffer with constant and vigorous agitation, 3 min each wash
<b>21 Quench Label Probe 6-AP</b> 35 min	<ul style="list-style-type: none"> <li>■ Add AP Stop QT to slides, incubate at RT in the dark for 30 min</li> <li>■ Wash slides 2 times in 1X PBS, 1 min each wash</li> <li>■ Rinse slides in Wash Buffer for 1 min</li> </ul>
<b>22 Label Probe 1-AP Hybridization</b> 20 min	<ul style="list-style-type: none"> <li>■ Prepare 1:1000 working Label Probe 1-AP solution in prewarmed Label Probe Diluent QF</li> <li>■ Add working Label Probe 1-AP solution to slides, incubate at 40 °C for 15 min</li> </ul>
<b>23 Wash Slides</b> 12 min	Wash slides 3 times in Wash Buffer, 3 min each wash
<b>24 Apply Fast Red Substrate</b> 45 min	<ul style="list-style-type: none"> <li>■ Add AP-Enhancer to slides and incubate at RT for 5 min</li> <li>■ Prepare Fast Red Substrate (1 Fast Red Tablet/5 mL Naphthol Buffer)</li> <li>■ Decant AP-Enhancer and add Fast Red Substrate to slides, incubate at 40 °C for 30 min</li> <li>■ Wash slides in 1X PBS for 1 min</li> </ul>

Step	Task
<b>25 Counterstain</b> <b>25 min</b>	<ul style="list-style-type: none"> <li>■ Incubate slides in Gill's Hematoxylin I stain at RT for 5-10 sec</li> <li>■ Wash slides 3 times in ddH<sub>2</sub>O, 1 min each wash</li> <li>■ Incubate in 0.01% ammonium hydroxide at RT for 10 sec</li> <li>■ Wash slides in ddH<sub>2</sub>O for 1 min</li> <li>■ Optional – Incubate slides in DAPI at RT for 1 min, wash slides in ddH<sub>2</sub>O for 1 min</li> <li>■ Let slides completely air dry at RT (~20 min)</li> </ul>
<b>26 Mount and Image</b> <b>40 min</b>	<div style="display: flex; justify-content: space-between;"> <div style="width: 60%;"> <p><b>DAKO Ultramount:</b></p> <p><b>For no coverslipping</b></p> <ul style="list-style-type: none"> <li>■ Add Ultramount to tissue sections</li> <li>■ Place slides in a 70 °C oven/incubator for 10-30 min</li> <li>■ Observe under brightfield or fluorescence microscope</li> <li>■ Store slides at RT</li> </ul> <p><b>For post mounting with coverslip</b></p> <ul style="list-style-type: none"> <li>■ Follow the no coverslip procedure</li> <li>■ Allow the slides to come to RT</li> <li>■ Apply Histomount directly on top of the dried Ultramount</li> <li>■ Place coverslip</li> <li>■ Air dry at RT for 15 min</li> <li>■ Observe under brightfield or fluorescence microscope</li> <li>■ Store slides at RT</li> </ul> </div> <div style="width: 35%;"> <p><b>Innovex Advantage Mounting Media:</b></p> <ul style="list-style-type: none"> <li>■ Add Advantage Mounting Media to cover glass and invert tissue slide to cover</li> <li>■ Flip over, allow slides to dry</li> <li>■ Seal all four edges with nail polish</li> <li>■ Observe under brightfield or fluorescence microscope</li> <li>■ Store slides at RT</li> </ul> </div> </div>

[www.affymetrix.com/panomics](http://www.affymetrix.com/panomics)

For research use only. Not for use in diagnostic procedures.