

# Quick Reference Guide

## QuantiGene ViewRNA miRNA ISH Cell Assay

### 24-Well Format for Adherent Cells

**!** **IMPORTANT:** If running QuantiGene ViewRNA miRNA ISH Cell Assay for the first time, please refer to the appropriate *QuantiGene ViewRNA miRNA ISH Cell Assay User Manual*, based on your sample type and processing format, to review assay guidelines and detailed procedures.

### Part 1: Sample Preparation

Step	Task
<b>Step 1. Coating Glass Coverslips with Poly-L-Lysine</b> 1 hr	<ul style="list-style-type: none"><li>■ Prepare 25 mL of 0.01% poly-L-lysine solution</li><li>■ Ensure availability of 100% ethanol, sterile 1X PBS, 24-well tissue culture plate, 12-mm glass coverslips, fine-tipped forceps</li><li>■ Sterilize coverslips in 24-well plate with 100% ethanol for 15 min at RT</li><li>■ Remove ethanol and air dry coverslips completely for 15 min at RT</li><li>■ Coat coverslips with diluted poly-L-lysine for 15 min at RT</li><li>■ Wash coverslips 3 times in 1X PBS</li><li>■ Air dry coated coverslips O/N at RT in laminar flow hood</li></ul>
<b>Step 2. Culturing Adherent Cells on Poly-L-Lysine Coated Coverslips</b> 20 min	<ul style="list-style-type: none"><li>■ Seed cells on coated coverslips</li><li>■ Culture cells O/N to 70-90% density</li></ul>

### Part 2: Pretreatment and Target Probe(s) Hybridization

Step	Task
<b>Step 3. Fix Samples</b> 1 hr 15 min	<ul style="list-style-type: none"><li>■ Prepare 1L of 1X PBS, 11 mL of 4% formaldehyde in 1X PBS</li><li>■ Wash cells 2 times with 1X PBS</li><li>■ Fix cells in 4% formaldehyde for 1 hr at RT (in a fume hood)</li><li>■ Remove formaldehyde and wash 3 times with 1X PBS</li></ul>
<b>Step 4. Prepare Buffers and Reagents during Sample Fixation</b> 10 min	<ul style="list-style-type: none"><li>■ Set a dry incubator to <math>40 \pm 1</math> °C</li><li>■ Equilibrate EDC to RT</li><li>■ Pre-warm Probe Set Diluent QF to 40 °C</li><li>■ Thaw Probe Set(s) and place on ice</li><li>■ Place Protease QS on ice</li><li>■ Prepare 800 mL of Wash Buffer, 30 mL Storage Buffer</li><li>■ Ensure availability of 50%, 70%, and 100% ethanol if dehydrating cells for storage or rehydrating cells for the assay</li><li>■ Ensure availability of Cross-Linking Buffer QM, Detergent Solution QC</li></ul>
<b>Step 5. Cross-Link with EDC</b> 1 hr 30 min	<ul style="list-style-type: none"><li>■ Equilibrate cells with Cross-Linking Buffer QM 2 times, 10 min each at RT with rocking</li><li>■ Prepare 0.16M EDC just before use</li><li>■ Cross-link target with 0.16 M EDC for 1 hr at RT, rocking</li><li>■ Wash cells 3 times with 1X PBS</li></ul>
<b>Step 6. Optional. Dehydrate/Rehydrate</b> 10-20 min	<ul style="list-style-type: none"><li>■ Dehydrate cells if storing sample: 50% EtOH 2 min -&gt; 70% EtOH, 2 min -&gt; 100% EtOH, 2 min -&gt; Fresh 100% EtOH and store at -20 °C</li><li>■ Rehydrate cells for assay if starting with dehydrated sample: 100% EtOH, 2 min -&gt; 70% EtOH, 2 min -&gt; 50% EtOH, 2 min -&gt; 1X PBS, 10 min</li></ul>

Step	Task
<b>Step 7. Permeabilize with Detergent and Digest with Protease</b> 40 min	<ul style="list-style-type: none"> <li>■ Permeabilize cells with Detergent Solution QC for 10 min at RT</li> <li>■ Wash cells 2 times with 1X PBS</li> <li>■ Prepare Working Protease Solution</li> <li>■ Treat cells with Working Protease Solution for 10 min at RT</li> <li>■ Wash cells 3X with 1X PBS</li> </ul>
<b>Step 8. Hybridize Target Probe(s)</b> 3 hrs	<ul style="list-style-type: none"> <li>■ Prepare Working Probe Set(s) Solutions</li> <li>■ Hybridize cells with Working Probe Solutions for 3 hrs at 40 °C</li> </ul>
<b>Step 9. Wash Samples</b> 10 min	<ul style="list-style-type: none"> <li>■ Wash cells 3 times in Wash Buffer, soaking 3 min each</li> </ul>
<b>Step 10. Stop Point</b> 1 min	<ul style="list-style-type: none"> <li>■ Store samples in Storage Buffer, covered with lid and seal with parafilm at 4 °C. Do not exceed 24 hrs</li> </ul>

### Part 3. Amplification and Detection

Step	Task
<b>Step 11. Prepare Additional Buffers and Reagents</b> 10 min	<ul style="list-style-type: none"> <li>■ Prewarm PreAmplifier Mix QM, Amplifier Mix QM, and Label Probe Diluent QF to 40 °C</li> <li>■ Equilibrate the following reagents to RT: <ul style="list-style-type: none"> <li>□ AP Enhancer, Fast Red Tablets, Naphthol Buffer</li> <li>□ Cell samples</li> </ul> </li> <li>■ Thaw, mix and place on ice 100X DAPI and Label Probe 4-Alexa 488</li> <li>■ Prepare 11 mL 1X DAPI and store on ice protected from light</li> <li>■ Place Label Probe 1-AP on ice</li> <li>■ Ensure availability of Wash Buffer and 1X PBS</li> <li>■ Ensure availability of Prolong® Gold Antifade Reagent</li> </ul>
<b>Step 12. Wash Storage Buffer off Samples</b> 5 min	<ul style="list-style-type: none"> <li>■ Wash cells 2X with Wash Buffer</li> </ul>
<b>Step 13. Hybridize with PreAmplifier</b> 1 hr	<ul style="list-style-type: none"> <li>■ Prepare Working PreAmplifier Mix Solution</li> <li>■ Hybridize cells with Working PreAmplifier Mix Solution for 1 hr at 40 °C</li> </ul>
<b>Step 14. Wash Samples</b> 10 min	<ul style="list-style-type: none"> <li>■ Wash cells 3X with Wash Buffer, soaking 3 min each</li> </ul>
<b>Step 15. Hybridize with Amplifier</b> 1 hr	<ul style="list-style-type: none"> <li>■ Prepare Working Amplifier Mix Solution</li> <li>■ Hybridize cells with Working Amplifier Mix Solution for 1 hr at 40 °C</li> </ul>
<b>Step 16. Wash Samples</b> 10 min	<ul style="list-style-type: none"> <li>■ Wash cells 3X with Wash Buffer, soaking 3 min each</li> </ul>
<b>Step 17. Hybridize with Label Probe(s)</b> 1 hr	<ul style="list-style-type: none"> <li>■ Prepare Label Probe(s) Working Solution</li> <li>■ Hybridize cells with Label Probe(s) Working Solution for 1 hr at 40 °C</li> </ul>
<b>Step 18. Wash Samples</b> 20 min	<ul style="list-style-type: none"> <li>■ Wash cells 3 times in Wash Buffer, soaking 3 min each for the first 2 washes and 10 min for the final wash</li> </ul>
<b>Step 19. Fast Red Development</b> 1 hr 20 min	<ul style="list-style-type: none"> <li>■ Incubate cells with AP-Enhancer for 5 min at RT; prepare Fast Red Substrate</li> <li>■ Remove AP-Enhancer, add Fast Red Substrate and incubate at 40 °C for 45 min</li> <li>■ Wash cells 2 times with 1X PBS</li> <li>■ Fix cells with 4% formaldehyde for 10 min at RT (in a fume hood)</li> <li>■ Wash cells 3 times with 1X PBS</li> </ul>

Step	Task
<b>Step 20. DAPI Counter Stain</b> <b>5 min</b>	<ul style="list-style-type: none"> <li>■ Stain cells with 1X DAPI for 5 min at RT</li> <li>■ Wash cells once with 1X PBS and replace with fresh 1X PBS</li> </ul>
<b>Step 21. Mount Samples on Glass Slides</b> <b>15 min</b>	<ul style="list-style-type: none"> <li>■ Mount glass coverslips to microscope glass slide</li> <li>■ Cure samples overnight at RT (protected from light)</li> <li>■ View under microscope                             <ul style="list-style-type: none"> <li>□ DAPI for nuclei: EX 387/11 nm, EM 447/60 nm (recommended Semrock P/N DAPI-11060B)</li> <li>□ ALEXA 488 for mRNA: EX 485/20 nm, EM 524/24 nm (recommended Semrock P/N FITC-2024B)</li> <li>□ Fast Red for miRNA: EX 530/40 nm, EM 590/40 nm (recommended Semrock P/N TRITC-B)</li> <li>□ ALEXA 750 for mRNA: EX 708/75 nm, EM 809/81 nm (recommended Semrock P/N Cy7-B)</li> </ul> </li> </ul>

[www.affymetrix.com/panomics](http://www.affymetrix.com/panomics)

For research use only. Not for use in diagnostic procedures.

P/N 19163 Rev. A 121004

© 2012 Affymetrix, Inc. All rights reserved. Affymetrix® and QuantiGene® are registered trademarks of Affymetrix, Inc. All other trademarks are the property of their respective owners.